



IDEAL CURES

Coating Troubleshooting Guide



Three Critical C's of Coating:

For effective troubleshooting of coating issues, delve into the Three C's. Identifying the precise root cause enables the production team or formulator to devise optimal solutions swiftly.

Core Tablet:

Physical attributes surface hardness and friability important considering stressful conditions of the coating process.

Coating Process:

Well optimised coating process parameters and ideal spray gun and equipment setup.

Coating Formulation:


Select an optimized formulation that aligns with core formulation requirements.



IDEAL CURES

INSTACOAT[®]
FILM COATING SYSTEMS

Tablet Breakage



Causes

- Poor core characteristics i.e. poor shape, too soft, brittle in nature
- Pan speed too high
- Inappropriate baffle design
- Pan load too low

Remedies

- Improve core formulation
- Reduce pan speed
- Pan load too low
- Improve baffle design
- Optimize pan load

Orange Peel



Causes

- Coating suspension viscosity too high
- Atomizing air pressure too low
- Spray rate too high

Remedies

- Optimize coating formula use an INSTACOAT high solids low viscosity system
- Increase atomizing air pressure
- Reduce the coating suspension solids level
- Decrease spray rate

Logo Infilling




Causes

- High gun to bed distance
- High atomizing air pressure
- Gun-to-bed distance too high
- Drying air temperature too high
- High solids content in the coating suspension

Remedies

- Decrease gun to bed distance
- Decrease atomizing air pressure
- Optimize spray gun set up
- Reduce drying air temperature
- Reduce coating suspension solids level

Tablet Discoloration



Causes

- Migration/ interaction of core/ coating ingredients often facilitated by residual solvent, use of liquid plasticizer, moisture uptake on storage.
- Fading due to poor pigment stability.
- Core defects not masked by coating.

Remedies

- Identify/ replace core ingredients
- Reduce over wetting by optimizing spray rate, drying air temperature & flow rate, pan speed, gun set-up, coating suspension solid level
- Identify/ replace coating ingredients
- Select an INSTACOAT moisture protective coating system
- Select stable pigments
- Select a higher opacity coating system

Tablet Edge Erosion




Causes

- Poor core tablet design i.e. sharp edges
- Worn or damaged tablet tooling
- Low film strength
- Low solids coating suspension
- High pan speed
- Low coating spray rate
- Low pan load

Remedies

- Improve core tablet design
- Replace or repair tooling
- Use a high film strength INSTACOAT system
- Use a high solids INSTACOAT system
- Reduce pan speed
- Optimize spray rate
- Optimize pan load

Picking & Sticking



Causes

- Coating solution spray rate too high
- Tablets bed temperature too low
- Inlet air volume too low
- Gun to bed distance too short
- Pan speed too low
- Coating suspension solids too low

Remedies

- Optimise the spray rate
- Increase coating bed temperature
- Increase inlet air volume
- Optimise gun to bed distance
- Increase pan speed
- Use a high solids INSTACOAT system

Logo Bridging



Causes

- Poor logo or break-line design (too narrow, too deep, square edges)
- Excipients with poor adhesion
- Overwet process conditions:
 - High suspension flow rate
 - Low gun- to-bed distance
 - Low pan Load
- Low adhesion coating formulation
- Change in solvents system from organic to aqueous

Remedies

- Improve logo design. Consult tablet tool supplier
- Choose high adhesion materials e.g. MCC, lactose
- Increase drying air temperature & flow rate
- Optimize gun position
- Increase tablet pan load
- Use INSTACOAT high adhesion coating system

Tablet Twinning



Causes

- Tablet shape i.e. flat shape, caplets
- Spray rate too high
- Gun to bed distance too close
- Low atomizing air pressure i.e. too big droplet size
- Coating process too wet

Remedies

- Change shape/ modify design of the core tablets i.e. avoid any perfectly flat areas
- Reduce spray rate
- Increase the gun to bed distance
- Increase the atomizing air pressure
- Increase the product bed/ exhaust temperature and / or drying air volume

Surface Erosion



Causes

- Poor tablet surface hardness i.e. friable core
- Tablets are too hygroscopic
- Logo design and placement
- Low film strength
- Low solids coating suspension
- Spray rate too low/ too high
- Pan speed too high

Remedies

- Modify core formulation
- Avoid or minimize super disintegrates in the core
- Modify logo design or placement
- Use INSTACOAT high strength film coating
- Optimize spray rate
- Decrease pan speed

Scuffing



Causes

- Tablet shape
- High tablet bed temperature
- Poor quality of stainless steel
- Improper cleaning of the coating pan
- High TiO_2 concentration in the flim

Remedies

- Avoid deep convex shapes
- Reduce drying air temperature
- Increase spray rate
- Polishing, passivating of stainless steel may help
- Pre-coat the pan with spray suspension
- Ensure proper cleaning of the coating pan
- Optimize TiO_2 concentration

Film Cracking/ Core Expansion




Causes

- Core has high rate of thermal expansion
- Post-compression core relaxation
- Higher coating bed temperature i.e. core expansion

Remedies

- Use organic fillers e.g. MCC, lactose, starch; avoid inorganics e.g. carbonates, phosphates
- Extend interval between tableting and coating process
- Decrease coating bed temperature

Dark Spots In Coating



Causes

- Poor pigment dispersion

Remedies

- Improve coating suspension preparation method (faster stirring, longer time)
- Use an Instacoat coating system

Film Splitting / Peeling



Causes

- Coating performed at high tablet bed temperature
- Poorly plasticized coating formulation
- Poor film strength coating formulation

Remedies

- Reduce bed temperature
- Use an INSTACOAT fully optimised coating formulation

“All rights reserved. No Part of photographs of image or design may be reproduced or copied.”

@ January 2024, IDEAL CURES PVT.LTD.

IDEAL CURES PRODUCT RANGE



IMMEDIATE RELEASE & MOISTURE BARRIER COATING

INSTACOAT 4G

INSTACOAT EHP 250

INSTACOAT AQUA II

INSTACOAT AQUA III

INSTACOAT UNIVERSAL

INSTACOAT AQUA LUSTER

INSTACOAT QD

INSTACOAT T2F

INSTACOAT SMART

INSTACOAT SOL

INSTACOAT ALTIMATE

INSTACOAT P4

INSTACOAT EMB

INSTAMOISTSHIELD
AQUA II

INSTAMOISTSHIELD

INSTANUTE MB II



Delayed Release Coating

INSTACOAT EN SUPER IV

INSTACOAT EN SUPER II

INSTACOAT HPMC- P

INSTACOAT EEN

INSTACOAT EEN SF



Pharma Acrylic Polymers

ECOPOL L30 D 55 (Liquid)

ECOPOL L100 55 (Dry Powder)

ECOPOL L100 (Dry Powder)

ECOPOL S100 (Dry Powder)



Extended Cooling Boosters

ECOCOOL
Pellets / Granules

ECOCOOL MP

ECOCOOL HD

ECOCOOL TI

ECOCOOL DT

ECOCOOL DP



Espheres InstaSpheres

ESPHERES EM

ESPHERES SiO₂

ESPHERES TA

MANNITOL PELLETS

INSTASPHERES EM

INSTASPHERES SiO₂

INSTASPHERES TA



Head Office and R&D Center II:
IDEAL CURES PVT. LTD.

A/ 223- 229, Virwani Industrial Estate, W. E. Highway,
Goregoan (East), Mumbai-400 063, India.



+91 022 42688700



info@idealcures.com



www.idealcures.com



Ideal Cures

Manufacturing Facility : Vasai • Jammu • Khambhat • Sikkim

INSTACOAT[®]
FILM COATING SYSTEMS





IDEAL CURES